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***Olivier Archer***

*Directeur Ventes & Marketing*

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*Coordinateur filiale France*

*Director New Business Development EMEA*

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## *“Contributing to Society through Healthcare”*

We contribute to society by providing valued products and services in the healthcare market and by responding to the needs of patients and healthcare professionals.

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

	TCD-10023PK <sup>1</sup>	CENTURY <sup>2,3</sup>	CENTURY II <sup>4,5</sup>	MASTER <sup>6</sup>	DISCOVERY 1TO3 <sup>7</sup>	CENTURY JSV <sup>8</sup>	e-Ultimaster <sup>9,10</sup>	MASTER DAPT <sup>11</sup>
<b>Nombre de patients</b>	22	105	1123	500	60	70	37,000	4300
<b>Design</b>	Single arm pharmacokinetics	Single arm, first-in-man study	Randomised 1:1 vs Xience	Randomised 3:1 vs BMS in patients with STEMI	Single arm, patients with multivessel disease	Single arm, patients requiring 2.25 mm diameter stents	Single arm, all-comers, real-world use; investigation of reduced DAPT regimens on clinical outcomes	Investigator-initiated, randomised, short vs prolonged DAPT
<b>Résultats primaires</b>	Sirolimus concentration in peripheral blood samples 28 days after Ultimaster implantation	Late loss at 6 months	Freedom from TLF at 9 months	Safety at 1 month, efficacy at 6 months, safety and efficacy at 12 months	OFDI strut coverage at 3 months	MACE at 9 months	TLF at 1 year	NACE, MACCE, MCB at 11 months (12 months post-index PCI)
<b>Etape</b>	Published	Published and completed	Published and completed	Submitted for publication, completed	Published	Published, Follow-up ongoing	Follow-up on going	Enrolling

Assessment of 'safety' refers to the side-effect profile and clinical events experienced (eg stroke, myocardial infarction, unplanned revascularisation, stent thrombosis). BMS, bare-metal stent; DAPT, dual antiplatelet therapy; MACCE, major adverse cardiac and cerebral events; MACE, major adverse coronary events; MCB, major or clinically relevant non-major bleeding; NACE, net adverse clinical endpoints; OFDI, optical frequency domain imaging; PCI, percutaneous coronary intervention; STEMI, ST-segment elevation myocardial infarction; TLF, target lesion failure.

1. Stojkovic S et al. *Fundam Clin Pharmacol* 2014;29:95–105; 2. Barbato E et al. *EuroIntervention* 2015;11:541–8; 3. Beleslin B. Presented at EuroPCR 2017, abstract OP0028; 4. Saito S et al. *Eur Heart J* 2014;35:2021–31; 5. Data on file at Terumo Corporation; 6. Valdes-Chavarrri M. Presented at PCR 2016; 7. Chevalier B et al. *Circ Cardiovasc Interv* 2017;10:pii:e004801;doi:10.1161/CIRCINTERVENTIONS.116.004801; 8. Saito S et al. *Cardiovasc Interv Ther* 2018;doi: 10.1007/s12928-018-0511-3; 9. e-ULTIMASTER trial. Available at: <https://clinicaltrials.gov/ct2/show/NCT02188355?term=e-ULTIMASTER&rank=1> (accessed April 2018); 10. Polad J. Presented at PCR2017, abstract POS0688; 11. MASTER DAPT trial. Available at <https://clinicaltrials.gov/ct2/show/NCT03023020?term=master+dapt&rank=1> (accessed April 2018).

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\*Compared with Ultimaster. †Passing through the second curve of a tortuous vessel model.

PCI, percutaneous coronary intervention.

Pushability is measured in terms of a force transmission through a vessel model. Kink resistance of the proximal shaft is assessed. Trackability is measured as the passing resistance through second curve of a model. Transmission force is assessed as a force transmission through a vessel model with unfixed guiding catheter. Bench test ISCD-523-31-18 performed by, and on file at, Terumo Corporation..


*Carotid Artery Stent*  
***Roadsaver***<sup>®</sup>

*Stent carotidien*

stent *Carotid Artery Stent* **Roadsaver**<sup>®</sup>

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